





passive engagement to reduce friction, improving the sliding mechanics during leveling and alignment stages, while provides effective torque and rotation control throughout the entire treatment, in its active mode. ZIP Lp is equipped with NiTi spring clip, strong enough for delivering continuous and consistent force and ensuring a tooth movements control. Its innovative "sandwich type" locking system prevents the clip from falling or blocking. Its retention is granted by an anatomical bondable 80 gauge mesh base, that increases its shear bond resistance reducing bonding failures.

FEATURES

- Low profile design for an enhanced patient's comfort;
- Lorem ipsum dolor sit amet, consectetuer adipiscing elit, sed diam nonummy nibh euismod tincidunt ut laoreet dolore.
- Easy opening/closing operations with a simple explorer for chairside time saving;
- Full slot area coverage for an outstanding tooth movement control throughout the treatment;
- 80 gauge mesh base, anatomically shaped for higher bond strength.

DESCRIPTION	ROTH 0.018"	ROTH 0.022"	MBT 0.018"	MBT 0.022"
Single Kit Hook 3/4/5	SLB30 - 11CCDH	SLB30 - 11CCVH	SLB30 - 14CCDH	SLB30 - 14CCVH
Single Kit Hook 3	SLB30 - 11CCDH3	SLB30 - 11CCVH3	SLB30 - 14CCDH3	SLB30 - 14CCVH3
Single Kit no Hook	SLB30 - 11CCDW	SLB30 - 11CCVW	SLB30 - 14CCDW	SLB30 - 14CCVW

SIA Vertriebspartner





